

# E.N.T CATALOG



Whatsapp & Call: +92 3034515151

Medical staff Order online , While we offer Distributors  
and businesses bulk orders with custom branding




# Table of Contents




Otology



Rhinology



Sinus Surgery



Larynx & Tonsil



Microlaryngology





Maxillofacial

## Otology

- Ear Polypus Forceps
- Ear Syringes
- Micro Ear Forceps
- Micro Ear Scissors
- Micro Ear Knives
- Micro Ear Hooks
- Micro Ear Needles
- Micro Ear Curettes
- Micro Ear Speculas
- Mastoid Rongeurs
- Tuning Forks

**BLACK COATED**



Straight, smooth



**HMI-ENT-001  
HMI-ENT-006**

Straight, serrated



**HMI-ENT-002  
HMI-ENT-007**

Bent upwards, serrated



**HMI-ENT-003  
HMI-ENT-008**

Right, serrated





**HMI-ENT-004  
HMI-ENT-009**

Left, serrated



**HMI-ENT-005  
HMI-ENT-010**



Straight, smooth



**HMI-ENT-011  
HMI-ENT-015**

Upwards, serrated



**HMI-ENT-012  
HMI-ENT-016**

Right, serrated





**HMI-ENT-013  
HMI-ENT-017**









Left, serrated



**HMI-ENT-014  
HMI-ENT-018**




BLACK  
COATED




003A

BLACK  
COATED

Wire bending forceps




MC GEE  
HMI-ENT-041-  
HMI-ENT-042




MC GEE  
HMI-ENT-039-  
HMI-ENT-040

Angled downward




0.8 x 3.5 mm  
HMI-ENT-039  
HMI-ENT-041

Angled downward




0.8 x 6.0 mm  
HMI-ENT-040  
HMI-ENT-042

Malleus nipper



HOUSE-DIETER  
HMI-ENT-047-  
HMI-ENT-050



HOUSE-DIETER  
HMI-ENT-043-  
HMI-ENT-046

Up cutting



HMI-ENT-043  
HMI-ENT-047

Down cutting



HMI-ENT-044  
HMI-ENT-048

Right cutting




HMI-ENT-045  
HMI-ENT-049

Left cutting




HMI-ENT-046  
HMI-ENT-050

Introducing Forceps, for teflon prosthesis





FULLER  
HMI-ENT-052



FULLER  
HMI-ENT-051

BLACK  
COATED



Straight

Bent upwards

Right



Left

HMI-ENT-053  
HMI-ENT-057

HMI-ENT-054  
HMI-ENT-058

HMI-ENT-055  
HMI-ENT-059

HMI-ENT-056  
HMI-ENT-060



Straight

Bent upwards

Right



Left

HMI-ENT-061  
HMI-ENT-065

HMI-ENT-062  
HMI-ENT-066

HMI-ENT-063  
HMI-ENT-067

HMI-ENT-064  
HMI-ENT-068



Straight

Bent upwards

Right

Left

HMI-ENT-069  
HMI-ENT-073

HMI-ENT-070  
HMI-ENT-074

HMI-ENT-071  
HMI-ENT-075

HMI-ENT-072  
HMI-ENT-076

Sickle knife



SCHUKNECHT  
HMI-ENT-077  
14.5 cm/5¾"

Roller knives

Vertical



SCHUKNECHT  
HMI-ENT-078  
15.5 cm/6"

Horizontal



SCHUKNECHT  
HMI-ENT-079  
15.5 cm/6"

Micro knife



ROSEN  
15 cm/6"

Diameter

Ø 1.0 mm  
Ø 1.4 mm  
Ø 1.5 mm  
Ø 1.6 mm  
Ø 2.0 mm  
Ø 2.3 mm  
Ø 2.5 mm  
Ø 2.6 mm  
Ø 2.8 mm  
Ø 3.0 mm

HMI-ENT-083 HMI-ENT-090  
HMI-ENT-084  
HMI-ENT-091  
HMI-ENT-085 HMI-ENT-092  
HMI-ENT-080 HMI-ENT-086 HMI-ENT-093  
HMI-ENT-081 HMI-ENT-094  
HMI-ENT-082 HMI-ENT-095  
HMI-ENT-087 HMI-ENT-096  
HMI-ENT-088 HMI-ENT-097  
HMI-ENT-089 HMI-ENT-098

45°

90°

Sickle knife



WULLSTEIN  
HMI-ENT-099  
15.5 cm/6"

Mucosa knife



TABB  
HMI-ENT-100  
15.5 cm/6"

Flap knife



PLESTER  
HMI-ENT-101  
15 cm/6"

Ear knives




SEXTON  
HMI-ENT-103  
17 cm/6¾"

Bayonet






SEXTON-GERZOG  
HMI-ENT-104  
17 cm/6¾"

Tympanum needles



Micro ear needles





Micro hooks, Picks

|        | Straight    | 25° Angled  |
|--------|-------------|-------------|
| 0.3 mm | HMI-ENT-167 | HMI-ENT-171 |
| 0.5 mm | HMI-ENT-168 | HMI-ENT-172 |
| 0.6 mm | HMI-ENT-169 | HMI-ENT-173 |
| 1.0 mm | HMI-ENT-170 | HMI-ENT-174 |

SHAMBAUGH  
16 cm/6¼"

45° Angled 90° Angled

|        |             |             |
|--------|-------------|-------------|
| 0.3 mm | HMI-ENT-175 | HMI-ENT-183 |
| 0.5 mm | HMI-ENT-176 | HMI-ENT-184 |
| 0.6 mm | HMI-ENT-177 | HMI-ENT-185 |
| 1.0 mm | HMI-ENT-178 | HMI-ENT-186 |
| 1.2 mm | HMI-ENT-179 | HMI-ENT-187 |
| 1.5 mm | HMI-ENT-180 | HMI-ENT-188 |
| 2.0 mm | HMI-ENT-181 | HMI-ENT-189 |
| 2.5 mm | HMI-ENT-182 | HMI-ENT-190 |

Micro otology hooks

|        | 25° Upwards | 25° Downwards |
|--------|-------------|---------------|
| 0.3 mm | HMI-ENT-191 | HMI-ENT-194   |
| 0.6 mm | HMI-ENT-192 | HMI-ENT-195   |
| 1.0 mm | HMI-ENT-193 | HMI-ENT-196   |

45° Upwards 45° Downwards

|        |             |             |
|--------|-------------|-------------|
| 0.3 mm | HMI-ENT-197 | HMI-ENT-200 |
| 0.6 mm | HMI-ENT-198 | HMI-ENT-201 |
| 1.0 mm | HMI-ENT-199 | HMI-ENT-202 |

90° Upwards 90° Downwards

|        |             |             |
|--------|-------------|-------------|
| 0.3 mm | HMI-ENT-203 | HMI-ENT-206 |
| 0.6 mm | HMI-ENT-204 | HMI-ENT-207 |
| 1.0 mm | HMI-ENT-205 | HMI-ENT-208 |

15.5 cm/6¼"

Micro dissectors

|                            | angled to right | angled to left |
|----------------------------|-----------------|----------------|
| HMI-ENT-209<br>15 cm/6"    |                 |                |
| HMI-ENT-210<br>15.5 cm/6¼" |                 |                |
| HMI-ENT-211<br>15 cm/6"    |                 |                |
| HMI-ENT-212<br>15.5 cm/6¼" |                 |                |

FISCH  
15.5 cm/6¼"

HOUSE  
15.5 cm/6¼"

Micro measuring rods

|                         |                         |
|-------------------------|-------------------------|
| HMI-ENT-213<br>Ø 3.5 mm | HMI-ENT-215<br>Ø 4.0 mm |
| HMI-ENT-214<br>Ø 4.5 mm | HMI-ENT-216<br>Ø 5.0 mm |


Foreign body levers

MOLDENHAUER  
HMI-ENT-217  
16 cm/6¼"

ZAUFAL  
HMI-ENT-218  
15 cm/6"


Periosteal elevators

6 mm



LEMPERT  
HMI-ENT-219  
20.5 cm/8"

6 mm



LEMPERT  
HMI-ENT-220  
18 cm/7"

Ear scoop  
sharp



LEMPERT  
21.5 cm/8½"

HMI-ENT-221  
Fig. 1

HMI-ENT-222  
Fig. 2

HMI-ENT-223  
Fig. 3

HMI-ENT-224  
Fig. 4

HMI-ENT-225  
Fig. 5

Foreign body instrument

Micro ear periosteal elevator



ROSEN  
HMI-ENT-226  
15 cm/6"



WEBER-LOCH  
HMI-ENT-227  
17 cm/6"



SHAPLEIGH  
15 cm/6"

HMI-ENT-228  
Fig. 1

HMI-ENT-229  
Fig. 2

Micro curettes

Ø 1.0 mm



HMI-ENT-230



Ø 1.2 mm

Ø 1.5 mm



HMI-ENT-231



Ø 1.8 mm

Ø 1.0 mm



HMI-ENT-232



Ø 1.2 mm




HOUSE  
14.5 cm/5¾"

Ear curretes




|       |  |             |             |             |             |             |             |             |
|-------|--|-------------|-------------|-------------|-------------|-------------|-------------|-------------|
| sharp |  | HMI-ENT-233 | HMI-ENT-234 | HMI-ENT-235 | HMI-ENT-236 | HMI-ENT-237 | HMI-ENT-238 | HMI-ENT-239 |
|       |  | 1.8 mm      | 2.1 mm      | 2.6 mm      | 3.2 mm      | 3.7 mm      | 4.1 mm      | 4.4 mm      |
|       |  | Fig. 000    | Fig. 00     | Fig. 0      | Fig. 1      | Fig. 2      | Fig. 3      | Fig. 4      |
| blunt |  | HMI-ENT-240 | HMI-ENT-241 | HMI-ENT-242 | HMI-ENT-243 | HMI-ENT-244 | HMI-ENT-245 | HMI-ENT-246 |
|       |  | 1.8 mm      | 2.1 mm      | 2.6 mm      | 3.2 mm      | 3.7 mm      | 4.1 mm      | 4.4 mm      |
|       |  | Fig. 000    | Fig. 00     | Fig. 0      | Fig. 1      | Fig. 2      | Fig. 3      | Fig. 4      |
| sharp |  | HMI-ENT-247 | HMI-ENT-248 | HMI-ENT-249 | HMI-ENT-250 | HMI-ENT-251 | HMI-ENT-252 | HMI-ENT-253 |
|       |  | 1.8 mm      | 2.1 mm      | 2.6 mm      | 3.2 mm      | 3.7 mm      | 4.1 mm      | 4.4 mm      |
|       |  | Fig. 000    | Fig. 00     | Fig. 0      | Fig. 1      | Fig. 2      | Fig. 3      | Fig. 4      |
| blunt |  | HMI-ENT-254 | HMI-ENT-255 | HMI-ENT-256 | HMI-ENT-257 | HMI-ENT-258 | HMI-ENT-259 | HMI-ENT-260 |
|       |  | 1.8 mm      | 2.1 mm      | 2.6 mm      | 3.2 mm      | 3.7 mm      | 4.1 mm      | 4.4 mm      |
|       |  | Fig. 000    | Fig. 00     | Fig. 0      | Fig. 1      | Fig. 2      | Fig. 3      | Fig. 4      |


Ear Scoops



Ear scoop sharp



Ear loops





HMI-ENT-268  
Ø 3 mm  
Fig. 0



HMI-ENT-269  
Ø 4 mm  
Fig. 1



HMI-ENT-270  
Ø 5 mm  
Fig. 2



HMI-ENT-271  
Ø 6 mm  
Fig. 3



HMI-ENT-272  
Ø 7.5 mm  
Fig. 4



HARTMANN  
3.6 cm / 1 1/2"

Fig.

|          |     |             |
|----------|-----|-------------|
| Set of 3 | 0-3 | HMI-ENT-273 |
| Set of 3 | 2-4 | HMI-ENT-274 |
| Set of 4 | 1-4 | HMI-ENT-274 |

Chrome plated



HARTMANN

Fig.

|          |            |             |
|----------|------------|-------------|
| Set of 3 | 1, 2, 4    | HMI-ENT-284 |
| Set of 4 | 3, 5, 7, 9 | HMI-ENT-285 |
| Set of 4 | 2, 4, 6, 8 | HMI-ENT-286 |



HMI-ENT-275  
Ø 3 mm  
Fig. 1



HMI-ENT-280  
Ø 4 mm  
Fig. 2



HMI-ENT-276  
Ø 4.5 mm  
Fig. 3



HMI-ENT-281  
Ø 5.0 mm  
Fig. 4



HMI-ENT-277  
Ø 5.5 mm  
Fig. 5



HMI-ENT-282  
Ø 6.0 mm  
Fig. 6



HMI-ENT-278  
Ø 6.5 mm  
Fig. 7



HMI-ENT-283  
Ø 7.0 mm  
Fig. 8



HMI-ENT-279  
Ø 7.5 mm  
Fig. 9



HMI-ENT-287  
4.5 x 5 mm



HMI-ENT-288  
5.5 x 6 mm



HMI-ENT-289  
6.5 x 7 mm



HMI-ENT-290  
7.5 x 8 mm

Chrome plated, Oval



HARTMANN



HMI-ENT-293  
Ø 4 mm  
Fig. 1



HMI-ENT-294  
Ø 5 mm  
Fig. 2



HMI-ENT-295  
Ø 6 mm  
Fig. 3



HMI-ENT-296  
Ø 7.5 mm  
Fig. 4



TOYNBEE  
3.6 cm / 1 1/2"

Fig.

|           |       |             |
|-----------|-------|-------------|
| Set for 3 | 2 - 4 | HMI-ENT-291 |
| Set for 4 | 1 - 4 | HMI-ENT-292 |



HMI-ENT-297  
Ø 8.2 mm



HMI-ENT-298  
Ø 7.7 mm



HMI-ENT-299  
Ø 7.2 mm



HMI-ENT-300  
Ø 6.7 mm



HMI-ENT-301  
Ø 6.3 mm



HMI-ENT-302  
Ø 5.8 mm



HMI-ENT-303  
Ø 5.3 mm



HMI-ENT-304  
Ø 4.8 mm



HMI-ENT-305  
Ø 4.3 mm



HMI-ENT-306  
Ø 3.8 mm



HMI-ENT-307  
Ø 3.3 mm



HMI-ENT-308  
9.0 mm



HMI-ENT-309  
8.5 mm



HMI-ENT-310  
8.0 mm



HMI-ENT-311  
7.5 mm



HMI-ENT-312  
7.0 mm



HMI-ENT-313  
6.5 mm



HMI-ENT-314  
6.0 mm



HMI-ENT-315  
5.5 mm



HMI-ENT-316  
5.0 mm



HMI-ENT-317  
4.5 mm



HMI-ENT-318  
4.0 mm



HMI-ENT-319  
9.0 mm



HMI-ENT-320  
8.5 mm



HMI-ENT-321  
8.0 mm



HMI-ENT-322  
7.5 mm



HMI-ENT-323  
7.0 mm



HMI-ENT-324  
6.5 mm



HMI-ENT-325  
6.0 mm



HMI-ENT-326  
5.5 mm



HMI-ENT-327  
5.0 mm



HMI-ENT-328  
4.5 mm




HMI-ENT-329  
4.0 mm

Satin finish, round



FARRIOR

Satin finish, oval



FARRIOR

Satin finish, oval



FARRIOR

HMI-ENT-344  
Ear specula storage rack



**HMI-ENT-330**  
Ø 11.5 mm  
Fig. 7



**HMI-ENT-331**  
Ø 11.5 mm  
Fig. 6



**HMI-ENT-332**  
Ø 10.5 mm  
Fig. 5



**HMI-ENT-333**  
Ø 8.5 mm  
Fig. 4



**HMI-ENT-334**  
Ø 7.5 mm  
Fig. 3



**HMI-ENT-335**  
Ø 6.5 mm  
Fig. 2



**HMI-ENT-336**  
Ø 5.5 mm  
Fig. 1



Round black

FARRIOR  
3.6 cm/1½"

HMI-ENT-343  
Fig. 1-6  
Complete set

Oval black



ZOELLNER  
3.8 cm/1½"



**HMI-ENT-337**  
11.5 mm  
Fig. 6



**HMI-ENT-338**  
10.5 mm  
Fig. 5



**HMI-ENT-339**  
9.5 mm  
Fig. 4



**HMI-ENT-340**  
8.5 mm  
Fig. 3



**HMI-ENT-341**  
7.5 mm  
Fig. 2



**HMI-ENT-342**  
6.5 mm  
Fig. 1



FRENCH  
HMI-ENT-345  
C<sup>2</sup> 16



FRENCH  
HMI-ENT-346  
C<sup>1</sup> 32



HARTMANN  
HMI-ENT-347  
C 64



HARTMANN  
HMI-ENT-348  
C 128



HARTMANN  
HMI-ENT-349  
C<sup>1</sup> 256



with dampers

ALUMINIUM



HMI-ENT-355

Set of 5  
(C 128 Hz - C 2048 Hz)

HMI-ENT-350  
C128

HMI-ENT-351  
C256

HMI-ENT-352  
C512

HMI-ENT-353  
C1024

HMI-ENT-354  
C2048



HARTMANN  
HMI-ENT-356  
C 128

HARTMANN  
HMI-ENT-357  
C¹ 256

HARTMANN  
HMI-ENT-358  
C² 512

HARTMANN  
HMI-ENT-359  
C³ 1024

HARTMANN  
HMI-ENT-360  
C⁴ 2048

HARTMANN  
HMI-ENT-361  
C⁵ 4096



LUCAE  
HMI-ENT-362  
C128



LUCAE  
HMI-ENT-363  
C1 256



LUCAE  
HMI-ENT-364  
C2 512



LUCAE  
HMI-ENT-366  
C3 1024

HMI-ENT-365  
Socket only



LUCAE  
HMI-ENT-367  
C4 2048



LUCAE  
HMI-ENT-368  
C5 4096



HMI-ENT-370  
Socket only

LUCAE  
HMI-ENT-369  
C 128



with hammer, pinwheel,  
brush & wheel for Neurc



RYDEL-SEIFFERT  
HMI-ENT-371  
C 64/C 128

HMI-ENT-372 HMI-ENT-374

Damper white scale for 40-4224-10  
Damper black scale for 40-4224-10

HMI-ENT-373 HMI-ENT-375

Knurled Screw for 40-4224-10  
Weights, C 64 for 40-4224-10

HMI-ENT-376  
C 64/C 128

KOENIG  
HMI-ENT-377  
A 109

40-1234-90  
Bakelit base



HMI-ENT-378  
Socket only










HARTMANN  
HMI-ENT-380  
C 128



HMI-ENT-381  
d 440

HARTMANN  
HMI-ENT-379  
C' 128








GOLDSMITH-THUDICHUM  
6.5 cm/2½"

18.5 x 9.5 mm



HMI-ENT-532  
Fig. 1

21.0 x 10.5 mm



HMI-ENT-533  
Fig. 2

21.0 x 11.5 mm



HMI-ENT-534  
Fig. 3



THUDICHUM (MARK HOVELL)

8.5 x 6.0 mm



HMI-ENT-535  
70 mm

9.0 x 7.0 mm



HMI-ENT-536  
71 mm

12.5 x 8.0 mm




HMI-ENT-537  
72 mm

15.0 x 8.5 mm




HMI-ENT-538  
73 mm

16.5 x 9.0 mm




HMI-ENT-539  
74 mm

18.0 x 10.0 mm




HMI-ENT-540  
74 mm

19.5 x 10.6 mm



HMI-ENT-541  
75 mm



Serrated



**HARTMANN-WULLSTEIN**  
1.8 x 4.0 mm HMI-ENT-403  
1.8 x 8.0 mm HMI-ENT-404

Smooth



**HARTMANN-WULLSTEIN**  
1.8 x 8.0 mm HMI-ENT-405



**HARTMANN-WULLSTEIN**  
8.5 cm/3½"  
Ø 1.0 mm HMI-ENT-406  
Ø 2.0 mm HMI-ENT-407  
Ø 2.5 mm HMI-ENT-408  
Ø 3.0 mm HMI-ENT-409

Serrated



**HARTMANN**  
8.5 cm/3¼" HMI-ENT-410  
14 cm/5½" HMI-ENT-411  
16.5 cm/6½" HMI-ENT-412

1:2



**HARTMANN**  
8.5 cm/3½" HMI-ENT-413  
14 cm/5½" HMI-ENT-414

Ø 3 mm



**HARTMANN**  
HMI-ENT-415  
8.5 cm/3¼"

Ø 2 mm



**HARTMANN-HERZFELD**  
HMI-ENT-416  
8.5 cm/3¼"

Ø 2.5 mm



**HARTMANN-HERZFELD**  
HMI-ENT-417  
8.5 cm/3¼"

Ø 3 mm



**HARTMANN-HERZFELD**  
HMI-ENT-418  
8.5 cm/3¼"

through cutting  
Ø 2 mm




**HARTMANN-HOFFMANN**  
HMI-ENT-419  
8.5 cm/3¼"




**HARTMANN-NOYES**  
Serrated HMI-ENT-420  
1x2 teeth HMI-ENT-421



**STRÜMPEL**  
HMI-ENT-422  
8.5 cm/3¼"



Mastoid Rongeurs




HMI-ENT-428  
3.5 x 3.5 mm

HMI-ENT-429  
4 x 4 mm

HMI-ENT-430  
4.5 x 4.5 mm

HMI-ENT-431  
5 x 5 mm

HMI-ENT-432  
6 x 6 mm



Valve Attachment



HMI-ENT-433

Suction Nozzle



HMI-ENT-434

Protecting Disc



HMI-ENT-435

Metal Olive



HMI-ENT-436

Conical tip



HMI-ENT-437

Button ended tip



HMI-ENT-438

Threaded tip



HMI-ENT-439



HMI-ENT-440

Flexible tube



HMI-ENT-441

Complete, self filling  
in plastic pouch

HMI-ENT-442

Rubber gasket

HMI-ENT-443

Retention spring

HMI-ENT-444

Syringe, only  
10 ml

## HMI-ENT-445

Complete ear syringe  
with shield & 2 tips

### Irrigation adapters



Shield

All metal syringe only Luer-Lock



### SCHIMMELBUSCH

|        |             |
|--------|-------------|
| 50 ml  | HMI-ENT-446 |
| 75 ml  | HMI-ENT-447 |
| 100 ml | HMI-ENT-448 |
| 150 ml | HMI-ENT-449 |
| 200 ml | HMI-ENT-450 |

### SCHIMMELBUSCH

|        |             |
|--------|-------------|
| 50 ml  | HMI-ENT-500 |
| 75 ml  | HMI-ENT-501 |
| 100 ml | HMI-ENT-502 |
| 150 ml | HMI-ENT-503 |
| 200 ml | HMI-ENT-504 |

### Accessories for SCHIMMELBUSCH



HMI-ENT-505  
Conical

HMI-ENT-506  
Olive

HMI-ENT-507  
Shield

HMI-ENT-508  
Olive large

HMI-ENT-509  
Conical small

HMI-ENT-510  
Conical large


HMI-ENT-511  
Olive bent

## Rhinology



- Adenoid Curettes
- Cartilage Crushers
- Chisels
- Elevators
- Osteotomes
- Gouges
- Nasal Cutting Forceps
- Nasal Dressing Forceps
- Nasal Knives
- Nasal Rasps
- Nasal Retractors
- Nasal Speculas
- Polypus Forceps
- Retractors
- Septum Forceps
- Skin Hooks



**HARTMANN**  
16.5 cm/6½"




**HARTMANN-HALLE**  
14 cm/5½"



**VIENNA**  
14 cm/5½"





**VIENNA**  
HMI-ENT-521  
14 cm/5½"





**TIECK-HALLE**  
HMI-ENT-522  
13 cm/5"



**KILLIAN**  
14 cm/5½"





GRUENWALD (JANSEN)



TROELTSCH  
HMI-ENT-544  
14.5 cm/5¾"



TROELTSCH  
HMI-ENT-545  
15.5 cm/6"



TROELTSCH  
15 cm/6" HMI-ENT-546  
18 cm/7" HMI-ENT-547



WESTMACOTT  
HMI-ENT-548  
20 cm/8"



LUBET-BARBON  
HMI-ENT-549  
18.5 cm/7¼"



HARTMANN  
HMI-ENT-550  
18.5 cm/7¼"

Nasal polypus forceps



**KNIGHT**  
HMI-ENT-551  
17 cm/6¾"



**GRAIG**  
HMI-ENT-552  
16 cm/6¼"




HMI-ENT-553



HMI-ENT-554



15.5 cm/6"




**LUC**  
20 cm/8"



HMI-ENT-555  
Fig. 1



HMI-ENT-556  
Fig. 2



**BRUENINGS**  
16.5 cm/6½"



HMI-ENT-557  
Fig. 1






HMI-ENT-558  
Fig. 2





HMI-ENT-559  
Fig. 3



HMI-ENT-560  
Fig. 4





## Nasal Cutting Forceps



HMI-ENT-578  
3 mm  
Fig. 1



HMI-ENT-579  
3.5 mm  
Fig. 2



HMI-ENT-580  
4.5 mm  
Fig. 3



HMI-ENT-581  
5.5 mm  
Fig. 4



HMI-ENT-582  
6 mm  
Fig. 5



HMI-ENT-583  
3 mm  
Fig. 1



HMI-ENT-584  
3.5 mm  
Fig. 2



HMI-ENT-585  
4.5 mm  
Fig. 3



HMI-ENT-586  
5.5 mm  
Fig. 4



HMI-ENT-587  
6 mm  
Fig. 5



HMI-ENT-588  
3 mm  
Fig. 1



HMI-ENT-589  
3.5 mm  
Fig. 2




HMI-ENT-590  
4.5 mm  
Fig. 3



HMI-ENT-591  
5.5 mm  
Fig. 4



HMI-ENT-592  
6 mm  
Fig. 5



WEIL-BLAKESLEY



HMI-ENT-593  
3 mm  
Fig. 1



HMI-ENT-594  
3.5 mm  
Fig. 2



HMI-ENT-595  
4.5 mm  
Fig. 3



HMI-ENT-596  
5.5 mm  
Fig. 4



HMI-ENT-597  
6 mm  
Fig. 5



HMI-ENT-598  
3 mm  
Fig. 1



HMI-ENT-599  
3.5 mm  
Fig. 2




HMI-ENT-600  
4.5 mm  
Fig. 3






HMI-ENT-601  
5.5 mm  
Fig. 4



HMI-ENT-602  
6 mm  
Fig. 5



BLAKESLEY





Ø 3 mm



½

KRAUSE  
HMI-ENT-628  
26 cm/10¼"



Ø 3 mm



½

KRAUSE-VOSS  
HMI-ENT-629  
26 cm/10¼"


Snare wire, coil of 10 m



½

Ø 0.3 mm HMI-ENT-630  
Ø 0.4 mm HMI-ENT-631

Snare Wire



½

HMI-ENT-632  
Ø 0.3 mm

Bundle with 12  
pieces

Through-cutting

TC Jaws




1/2  
MIDDLETON-JANSEN  
20.5 cm/8 1/4"  
4 x 20 mm HMI-ENT-633  
5 x 24 mm HMI-ENT-634



1/2  
MIDDLETON-JANSEN  
HMI-ENT-635  
20 cm/8"



1/2  
HMI-ENT-636  
20 cm/8"



RUBIN  
HMI-ENT-637  
20 cm/8"



RUBIN  
HMI-ENT-638  
18.5 cm/7 1/4"



RUBIN  
HMI-ENT-639  
18.5 cm/7 1/4"

Protection Cap



RUBIN  
HMI-ENT-640

One toothed cutting edge



CAPLAN  
HMI-ENT-641  
20 cm/8"


Longitudinal Jaw




HMI-ENT-642

TC Jaws


Morselizer Jaw




HMI-ENT-643




GUNTER  
18 cm/7"




TC Edge





CAPLAN  
HMI-ENT-644  
21 cm/8 1/4"




TC Edge



GORNEY  
HMI-ENT-645  
20.5 cm/8"





MCKENTY  
15 cm/6"




HALLE  
HMI-ENT-666  
16 cm/6¼"



COTTLE  
HMI-ENT-667  
19 cm/7½"




JOSEPH  
HMI-ENT-668  
16.5 cm/6½"



FREER  
19 cm/7½"


Cleft Palate Periosteal Rasapatory




FREER ½  
HMI-ENT-671  
18 cm/7"

FREER ½  
HMI-ENT-672  
18 cm/7"

FREER  
HMI-ENT-673  
21 cm/8¼"




KRENKEL  
24.5 cm/9½"




MCKENTY  
HMI-ENT-676  
20.5 cm/8"

BALLENGER-HAJEK  
HMI-ENT-677  
19 cm/7½"




**HOWARD**  
HMI-ENT-678  
21.5 cm/8½"




**COTTLE**  
HMI-ENT-679  
23.5 cm/9¼"



**KILLIAN**  
HMI-ENT-680  
18 cm/7"



**KILLIAN**  
HMI-ENT-681  
23.5 cm/9¼"




**PENNINGTON**  
HMI-ENT-682  
21 cm/8"




**ROGER**  
HMI-ENT-683  
19 cm/7½"



**JOSEPH-MASING**  
HMI-ENT-684  
21.5 cm/8½"



**COTTLE**  
HMI-ENT-685  
23 cm/9"



**COTTLE**  
HMI-ENT-686  
23 cm/9"



**JOSEPH**  
HMI-ENT-687  
18 cm/7"



FREER  
HMI-ENT-688  
16.5 cm/6½"



FREER  
HMI-ENT-689  
16 cm/6¼"




KILLIAN-CLAUS  
HMI-ENT-690  
16 cm/6¼"




KILLIAN  
HMI-ENT-691  
16 cm/6¼"



BALLENGER  
15 cm/6"




HMI-ENT-692  
6 mm




HMI-ENT-694  
6 mm



HMI-ENT-693  
8 mm



HMI-ENT-695  
6 mm



SHEEHAN  
16 cm/6¼"



RUBIN  
18.5 cm/7¼"



COTTLE  
18.5 cm/7¼"




COTTLE  
HMI-ENT-725  
18.5 cm/7¼"

|      |             |       |             |
|------|-------------|-------|-------------|
| 2 mm | HMI-ENT-696 | 9 mm  | HMI-ENT-703 |
| 3 mm | HMI-ENT-697 | 10 mm | HMI-ENT-704 |
| 4 mm | HMI-ENT-698 | 11 mm | HMI-ENT-705 |
| 5 mm | HMI-ENT-699 | 12 mm | HMI-ENT-706 |
| 6 mm | HMI-ENT-700 | 13 mm | HMI-ENT-707 |
| 7 mm | HMI-ENT-701 | 14 mm | HMI-ENT-708 |
| 8 mm | HMI-ENT-702 | 15 mm | HMI-ENT-709 |

|       |             |
|-------|-------------|
| 6 mm  | HMI-ENT-710 |
| 8 mm  | HMI-ENT-711 |
| 10 mm | HMI-ENT-712 |
| 11 mm | HMI-ENT-713 |
| 12 mm | HMI-ENT-714 |
| 14 mm | HMI-ENT-715 |
| 16 mm | HMI-ENT-716 |

|       |             |
|-------|-------------|
| 2 mm  | HMI-ENT-717 |
| 4 mm  | HMI-ENT-718 |
| 6 mm  | HMI-ENT-719 |
| 7 mm  | HMI-ENT-720 |
| 8 mm  | HMI-ENT-721 |
| 9 mm  | HMI-ENT-722 |
| 10 mm | HMI-ENT-723 |
| 12 mm | HMI-ENT-724 |



5 mm

Straight

Right

Left



HMI-ENT-749

HMI-ENT-750

HMI-ENT-751

NEIVERT-ANDERSON  
20.5 cm/8"

5 mm

Straight

Right

Left



HMI-ENT-752

HMI-ENT-753

HMI-ENT-754

SILVER  
18 cm/7"




PARKES  
HMI-ENT-755  
18 cm/7"


Nasal Saws



JOSEPH




HMI-ENT-756  
19 cm/7½"




HMI-ENT-757  
19 cm/7½"




HMI-ENT-758  
19 cm/7½"





HMI-ENT-759  
17 cm/6¾"




HMI-ENT-760  
17 cm/6¾"



MALTZ  
17 cm/6¾"



LEWIS  
19 cm/7½"




45 x 6 mm

GALLAHER  
17 cm/6¾"




HMI-ENT-765 HMI-ENT-766 HMI-ENT-767

34 x 8 mm




COTTLE  
HMI-ENT-768  
21 cm/8¼"



32 x 9 mm



AUFRICHT  
21 cm/8¼"




AUFRICHT (GLABELLA)  
20 cm/8"

Downward Upward  
4 mm HMI-ENT-771 HMI-ENT-773  
7 mm HMI-ENT-772 HMI-ENT-774


TC Jaws



AUFRICHT (GLABELLA)  
HMI-ENT-775  
21 cm/8¼"



JOSEPH  
HMI-ENT-776  
16 cm/6 1/4"



FOMON  
21.5 cm/8 1/2"  
Fine Coarse  
HMI-ENT-777 HMI-ENT-778



FOMON  
HMI-ENT-779  
20.5 cm/8"

Guide Support for osteotomes

Cartilage Crusher



COTTLE  
Smooth Bed HMI-ENT-780  
Serrated Bed HMI-ENT-781



TOLSDORF  
Smooth Bed HMI-ENT-782  
Serrated Bed HMI-ENT-783



COTTLE  
HMI-ENT-784  
7 mm




COTTLE  
HMI-ENT-785  
12 cm/4 3/4"

## TC Jaws




FOMON  
20 cm/8"

very coarse




1/2  
Teeth 12/14



HMI-ENT-786

coarse




3/4  
Teeth 16/20




HMI-ENT-787

medium




5/6  
Teeth 24/30




HMI-ENT-788

medium fine



7/8  
Teeth 36/44



HMI-ENT-789

very fine



9/10  
Teeth 52/62



HMI-ENT-790

## TC Jaws



PARKES  
21 cm/8 1/4"

very coarse



12

HMI-ENT-791  
Fig. 1

coarse



14

HMI-ENT-792  
Fig. 2



16

HMI-ENT-793  
Fig. 3



20

HMI-ENT-794  
Fig. 4

medium



24

HMI-ENT-795  
Fig. 5

medium



30

HMI-ENT-796  
Fig. 6

medium fine



36

HMI-ENT-797  
Fig. 7

fine



44

HMI-ENT-798  
Fig. 8

very fine




52


HMI-ENT-799  
Fig. 9




62

HMI-ENT-800  
Fig. 10





COTTLE-NEIVERT  
HMI-ENT-821  
19.5 cm/7 1/4"




HMI-ENT-822  
sharp

HMI-ENT-823  
blunt/sharp



HMI-ENT-824  
sharp/blunt

COTTLE  
19.5 cm/7 3/4"



COTTLE  
HMI-ENT-825  
14 cm/5 1/2"



NEIVERT  
HMI-ENT-826  
17 cm/6 3/4"

Saddle Hook



COTTLE  
HMI-ENT-827  
14 cm/5 1/2"



COTTLE  
HMI-ENT-828  
16 cm/6 1/4"

60 x 6.5 mm



AUFRICHT  
HMI-ENT-829  
12 cm/4 3/4"

4 mm



OGURA  
HMI-ENT-830  
11 cm/4 1/4"

1.7 mm



COTTLE  
HMI-ENT-831  
14 cm/5 1/2"

Columella Forceps

12 mm



COTTLE  
HMI-ENT-832  
11 cm/4 1/4"




COAKLEY  
17 cm/6 3/4"



HMI-ENT-833  
Fig. 1



HMI-ENT-834  
Fig. 2




HMI-ENT-835  
Fig. 3



HMI-ENT-836  
Fig. 4



HMI-ENT-837  
Fig. 5



HMI-ENT-838  
Fig. 6



GORNEY  
19 cm/7½"

3 mm HMI-ENT-839  
4 mm HMI-ENT-840



GUILLEN  
HMI-ENT-841  
19 cm/7½"

HMI-ENT-842  
Obturator  
for 41-1840-20



HMI-ENT-843  
Ø 3 mm

HMI-ENT-844  
Ø 4 mm

RITTER-HALLE  
19 cm/7½"



FAULKNER  
HMI-ENT-845  
21 cm/8¼"



BECKMANN  
½ 22 cm/8¾"

Fig.

|   |       |             |
|---|-------|-------------|
| 1 | 13 mm | HMI-ENT-846 |
| 2 | 15 mm | HMI-ENT-847 |
| 3 | 17 mm | HMI-ENT-848 |
| 4 | 19 mm | HMI-ENT-849 |
| 5 | 21 mm | HMI-ENT-850 |
| 6 | 23 mm | HMI-ENT-851 |



BECKMANN  
½ 21.5 cm/8½"

Fig.

|   |       |             |
|---|-------|-------------|
| 1 | 13 mm | HMI-ENT-852 |
| 2 | 15 mm | HMI-ENT-853 |
| 3 | 17 mm | HMI-ENT-854 |
| 4 | 19 mm | HMI-ENT-855 |
| 5 | 21 mm | HMI-ENT-856 |
| 6 | 23 mm | HMI-ENT-857 |



KRESSNER  
HMI-ENT-858  
22 cm/8¾"




COTTLE-WALSHAM  
16 cm/6¼"



HMI-ENT-859  
Straight



HMI-ENT-860  
Left




HMI-ENT-861  
Right





ASCH  
HMI-ENT-862  
22.5 cm/9"


## Sinus Surgery


- Antrum Cannula
- Antrum Curettes
- Antrum Punches
- Nasal Cutting Forceps
- Nasal Dressing Forceps
- Nasal Scissors
- Punches
- Sickle Knives
- Sinus Probes
- Sinus Seekers
- Spoon Forceps











## Larynx, Tonsil

- Injection Needles
- Laryngeal Forceps
- Laryngeal Mirrors
- Laryngeal Polypus Forceps
- Retractors
- Tampon Forceps
- Tongue Depressor
- Tonsil Hemostatic Forceps
- Tonsil Knives
- Tonsil Punches
- Tonsil Seizing Forceps
- Tracheal Instruments




Tonsil grasping forceps



WHITE  
HMI-ENT-936  
18 cm/7"




BLOHMKE  
HMI-ENT-937  
20 cm/8"




HMI-ENT-938

HMI-ENT-939


COLVER  
19 cm/7½"



SCHNIDT  
HMI-ENT-940



SCHNIDT  
HMI-ENT-942



SAWTELL  
HMI-ENT-941




SAWTELL  
HMI-ENT-943





HMI-ENT-944  
Fig. 1



HMI-ENT-945  
Fig. 2



TIVNEN-TYDINGS  
21 cm/8¼"





**BRUENINGS**  
22 cm/8 $\frac{1}{4}$ " HMI-ENT-957  
32 cm/12 $\frac{3}{4}$ " HMI-ENT-958  
41 cm/16 $\frac{1}{4}$ " HMI-ENT-959



**BRUENINGS**  
HMI-ENT-960  
18 cm/7"



**HUBER**  
HMI-ENT-961



**FRAENKEL**  
HMI-ENT-962



**SEIFFERT**  
HMI-ENT-963



**FRAENKEL**  
HMI-ENT-964



**SCHEINMANN**  
HMI-ENT-965



**SCHEINMANN**  
HMI-ENT-966



**BRUENINGS**  
Fig. 2 HMI-ENT-967  
Fig. 3 HMI-ENT-968



**BRUENINGS**  
Fig. 2 HMI-ENT-969  
Fig. 3 HMI-ENT-970




**LANGE**  
HMI-ENT-971



**KRAUSE**  
Fig. 1 HMI-ENT-972  
Fig. 2 HMI-ENT-973



**KRAUSE**  
HMI-ENT-974



**CORDES**  
Fig. 1 HMI-ENT-975  
Fig. 2 HMI-ENT-976  
Fig. 3 HMI-ENT-977  
Fig. 4 HMI-ENT-978



**SCHUMACHER**  
HMI-ENT-979



**CORDES**  
Fig. 1 HMI-ENT-980  
Fig. 2 HMI-ENT-981



**CORDES**  
HMI-ENT-982




**ROSENBERG**  
HMI-ENT-983

Plane, with fixed handle



18.5 cm/7 1/4"



|         |         |             |
|---------|---------|-------------|
| Fig. 00 | Ø 8 mm  | HMI-ENT-988 |
| Fig. 0  | Ø 10 mm | HMI-ENT-989 |
| Fig. 1  | Ø 12 mm | HMI-ENT-990 |
| Fig. 2  | Ø 14 mm | HMI-ENT-991 |
| Fig. 3  | Ø 16 mm | HMI-ENT-992 |
| Fig. 4  | Ø 18 mm | HMI-ENT-993 |
| Fig. 5  | Ø 20 mm | HMI-ENT-994 |
| Fig. 6  | Ø 22 mm | HMI-ENT-995 |
| Fig. 7  | Ø 24 mm | HMI-ENT-996 |
| Fig. 8  | Ø 26 mm | HMI-ENT-997 |
| Fig. 9  | Ø 28 mm | HMI-ENT-998 |
| Fig. 10 | Ø 30 mm | HMI-ENT-999 |


Moveable mirror




MICHEL  
27 cm/7 1/2"




HMI-ENT-984  
Fig. 1



HMI-ENT-985  
Fig. 2



HMI-ENT-986  
Fig. 3



HMI-ENT-987  
Fig. 4


Rhodium front coated

|         |         |              |
|---------|---------|--------------|
| Fig. 00 | Ø 8 mm  | HMI-ENT-1000 |
| Fig. 0  | Ø 10 mm | HMI-ENT-1001 |
| Fig. 1  | Ø 12 mm | HMI-ENT-1002 |
| Fig. 2  | Ø 14 mm | HMI-ENT-1003 |
| Fig. 3  | Ø 16 mm | HMI-ENT-1004 |
| Fig. 4  | Ø 18 mm | HMI-ENT-1005 |
| Fig. 5  | Ø 20 mm | HMI-ENT-1006 |
| Fig. 6  | Ø 22 mm | HMI-ENT-1007 |
| Fig. 7  | Ø 24 mm | HMI-ENT-1008 |
| Fig. 8  | Ø 26 mm | HMI-ENT-1009 |
| Fig. 9  | Ø 28 mm | HMI-ENT-1010 |
| Fig. 10 | Ø 30 mm | HMI-ENT-1011 |

Plane, without handle



13.5 cm/5 1/4"



|         |         |              |
|---------|---------|--------------|
| Fig. 00 | Ø 8 mm  | HMI-ENT-1012 |
| Fig. 0  | Ø 10 mm | HMI-ENT-1013 |
| Fig. 1  | Ø 12 mm | HMI-ENT-1014 |
| Fig. 2  | Ø 14 mm | HMI-ENT-1015 |
| Fig. 3  | Ø 16 mm | HMI-ENT-1016 |
| Fig. 4  | Ø 18 mm | HMI-ENT-1017 |
| Fig. 5  | Ø 20 mm | HMI-ENT-1018 |
| Fig. 6  | Ø 22 mm | HMI-ENT-1019 |
| Fig. 7  | Ø 24 mm | HMI-ENT-1020 |
| Fig. 8  | Ø 26 mm | HMI-ENT-1021 |
| Fig. 9  | Ø 28 mm | HMI-ENT-1022 |
| Fig. 10 | Ø 30 mm | HMI-ENT-1023 |

Rhodium front coated

|         |         |              |
|---------|---------|--------------|
| Fig. 00 | Ø 8 mm  | HMI-ENT-1024 |
| Fig. 0  | Ø 10 mm | HMI-ENT-1025 |
| Fig. 1  | Ø 12 mm | HMI-ENT-1026 |
| Fig. 2  | Ø 14 mm | HMI-ENT-1027 |
| Fig. 3  | Ø 16 mm | HMI-ENT-1028 |
| Fig. 4  | Ø 18 mm | HMI-ENT-1029 |
| Fig. 5  | Ø 20 mm | HMI-ENT-1030 |
| Fig. 6  | Ø 22 mm | HMI-ENT-1031 |
| Fig. 7  | Ø 24 mm | HMI-ENT-1032 |
| Fig. 8  | Ø 26 mm | HMI-ENT-1033 |
| Fig. 9  | Ø 28 mm | HMI-ENT-1034 |
| Fig. 10 | Ø 30 mm | HMI-ENT-1035 |




HMI-ENT-1036  
11 cm/4 1/4"



HMI-ENT-1037  
12.5 cm/5"



FRAENKEL  
HMI-ENT-1038  
18 cm/7"



½  
JURASZ  
23 cm/9"



HMI-ENT-1039  
Fig. 1



HMI-ENT-1040  
Fig. 2



HMI-ENT-1041  
Fig. 3




HMI-ENT-1042  
Fig. 4


13.5 mm



NAGER  
HMI-ENT-1043  
19 cm/17½"



Velum Hook



Arch palate hook

HOFER  
HMI-ENT-1044  
24 cm/9½"



REICHERT  
HMI-ENT-1045  
24 cm/9½"

3 mm



JACKSON  
HMI-ENT-1046  
22 cm/8¾"

3 mm




HMI-ENT-1047




## Microlaryngology


- Chest Supports
- Cotton Carriers
- Injection Cannulas
- Light Guide Carriers
- Microlaryngology Forceps
- Microlaryngology Instruments
- Microlaryngology Needle Holders
- Microlaryngology Scissors
- Operating Laryngoscopes
- Suction Tubes
- Teeth Protectors




HMI-ENT-1065  
17 cm/6¾"  
Adult X - Large



HMI-ENT-1066  
17 cm/6¾"  
Adult Large




HMI-ENT-1067  
17 cm/6¾"  
Adult Medium




HMI-ENT-1068  
18 cm/7"  
Adult Small


most commonly used model particularly difficult anatomical circumstances




HMI-ENT-1069  
18 cm/7"  
Adult Medium



HMI-ENT-1070  
22 cm/8¾"  
Adult Medium



HMI-ENT-1071  
15 cm/6"  
Children Small



HMI-ENT-1072  
13 cm/5"  
Infant Small

#### Fiber Optic Light Guide Clip

For proximal illumination

HMI-ENT-1073

#### Fiber Optic Light Guide Carriers

for distal illumination

|              | Working length | Use with                    |
|--------------|----------------|-----------------------------|
| HMI-ENT-1074 | 7.5 cm/3"      | HMI-ENT-1071 - HMI-ENT-1072 |
| HMI-ENT-1075 | 14 cm/5½"      | HMI-ENT-1065- HMI-ENT-1070  |

#### Suction Tubes



To remove vapor, for LASER treatment, Dia. Ø 3 mm


|              | Working length | Use with                    |
|--------------|----------------|-----------------------------|
| HMI-ENT-1076 | 9 cm/3½"       | HMI-ENT-1071 - HMI-ENT-1072 |
| HMI-ENT-1077 | 12.5 cm/5"     | HMI-ENT-1065- HMI-ENT-1070  |

#### Injection Cannula


for positive pressure assisted ventilation

|              | Working length | Use with                    |
|--------------|----------------|-----------------------------|
| HMI-ENT-1078 | 12.5 cm/5"     | HMI-ENT-1071 - HMI-ENT-1072 |







Spreadable




HMI-ENT-1082  
17 cm/6¾"  
Adult X - Large



HMI-ENT-1083  
17 cm/6¾"  
Adult Large




HMI-ENT-1084  
17 cm/6¾"  
Adult Medium




HMI-ENT-1085  
18 cm/7"  
Adult Small


most commonly used model particularly difficult anatomical circumstances




HMI-ENT-1086  
18 cm/7"  
Adult Medium



HMI-ENT-1087  
22 cm/8¾"  
Adult Medium



HMI-ENT-1088  
15 cm/6"  
Children Small



HMI-ENT-1089  
13 cm/5"  
Infant Small

#### Fiber Optic Light Guide Clip

For proximal illumination

HMI-ENT-1090

#### Fiber Optic Light Guide Carriers

for distal illumination

|              | Working length | Use with                   |
|--------------|----------------|----------------------------|
| HMI-ENT-1091 | 7.5 cm/3"      | HMI-ENT-1071 -HMI-ENT-1072 |
| HMI-ENT-1092 | 14 cm/5½"      | HMI-ENT-1065- HMI-ENT-1070 |

#### Suction Tubes

To remove vapor, for LASER treatment, Dia. Ø 3 mm

|              | Working length | Use with                   |
|--------------|----------------|----------------------------|
| HMI-ENT-1093 | 9 cm/3½"       | HMI-ENT-1071 -HMI-ENT-1072 |
| HMI-ENT-1094 | 12.5 cm/5"     | HMI-ENT-1065- HMI-ENT-1070 |


#### Injection Cannula

for positive pressure assisted ventilation


|              | Working length | Use with                   |
|--------------|----------------|----------------------------|
| HMI-ENT-1095 | 12.5 cm/5"     | HMI-ENT-1071 -HMI-ENT-1072 |



HMI-ENT-1096



HMI-ENT-1097



HMI-ENT-1098

#### KLEINSASSER Chest Support

With adjustment wheel, Ring Diameter Ø 9.5 mm

|              | Working length | for      |
|--------------|----------------|----------|
| HMI-ENT-1099 | 24 cm/9½"      | Children |
| HMI-ENT-1100 | 34 cm/13½"     | Adults   |



with adjustment wheel

#### Support Rod Only, With ring

|              | Working length | for      |
|--------------|----------------|----------|
| HMI-ENT-1101 | 24 cm/9½"      | Children |
| HMI-ENT-1102 | 34 cm/13½"     | Adults   |



#### Holder Only


|              | for           |
|--------------|---------------|
| HMI-ENT-1103 | Chest Support |






**HOLINGER Operating Laryngoscope**  
Complete with light guide

|              | Use for  | Handles | Working Length |
|--------------|----------|---------|----------------|
| HMI-ENT-1104 | Babies   | Small   | 12.5 cm/5"     |
| HMI-ENT-1105 | Children | Small   | 17 cm/6¾"      |
| HMI-ENT-1106 | Adult    | Large   | 18 cm/7"       |



**JACKSON Operating Laryngoscope**  
Complete with light guide

|              | Use for  | Handles | Working Length |
|--------------|----------|---------|----------------|
| HMI-ENT-1107 | Babies   | Small   | 12.5 cm/5"     |
| HMI-ENT-1108 | Children | Small   | 17 cm/6¾"      |
| HMI-ENT-1109 | Adult    | Large   | 18 cm/7"       |




MicroLaryngology Forceps - Finger Movable





| Ø 2 mm       | Grasping Forceps |  |
|--------------|------------------|--|
| HMI-ENT-1110 | Straight         |  |
| HMI-ENT-1111 | Upwards          |  |
| HMI-ENT-1112 | Curved right     |  |
| HMI-ENT-1113 | Curved left      |  |

| Ø 2 mm       | Scissors          |  |
|--------------|-------------------|--|
| HMI-ENT-1114 | Straight          |  |
| HMI-ENT-1115 | Upwards 45°       |  |
| HMI-ENT-1116 | Curved right      |  |
| HMI-ENT-1117 | Curved left       |  |
| HMI-ENT-1118 | Curved horizontal |  |

| Ø 4 mm       | Ø 2 mm       | Cup Shaped Forceps |  |
|--------------|--------------|--------------------|--|
| HMI-ENT-1119 | HMI-ENT-1120 | Straight           |  |
| HMI-ENT-1121 | HMI-ENT-1122 | Upwards            |  |
| HMI-ENT-1123 | HMI-ENT-1124 | Curved right       |  |
| HMI-ENT-1125 | HMI-ENT-1126 | Curved left        |  |



|              |                                |  |
|--------------|--------------------------------|--|
|              |                                |  |
| HMI-ENT-1128 | knife, oval, straight          |  |
| HMI-ENT-1129 | knife, oval, 45°               |  |
| HMI-ENT-1130 | sickel knife, curved           |  |
| HMI-ENT-1131 | sickel knife, straight         |  |
| HMI-ENT-1132 | knife, golf stick shaped       |  |
| HMI-ENT-1133 | knife, round, vertical cutting |  |
| HMI-ENT-1134 | hook, blunt                    |  |
| HMI-ENT-1135 | hook, sharp                    |  |
| HMI-ENT-1136 | hook, with probe end 90°       |  |
| HMI-ENT-1137 | needle, curved to right        |  |
| HMI-ENT-1138 | needle, curved to left         |  |
| HMI-ENT-1139 | elevator, with suction channel |  |
| HMI-ENT-1140 | knot pusher / tier             |  |



## Oral & Maxilo Facial

- Bone Contouring Forceps
- Cheek Retractor
- Chest Support
- Extra Oral Retractors
- Intra Oral Retractors
- Lip & Cheek Vestibulum Retractors
- Mandibular Retractors
- Mouthgags
- Multi Purpose Retractors
- Osteotomes
- Periosteal Elevators
- Ramus Retractors
- Soft Tissue Retractors



43-4224-04  
Fram only



DINGMANN  
HMI-ENT-1148  
Complete Set

HMI-ENT-1147  
Fig. 1

HMI-ENT-1146  
Fig. 2

HMI-ENT-1145  
Fig. 3




McIVOR  
HMI-ENT-1152  
Complete set

HMI-ENT-1151  
Fig. 1



HMI-ENT-1150  
Fig. 2

HMI-ENT-1149  
Fig. 3



MOLT  
HMI-ENT-1160  
11 cm/5 1/2"

MOLT  
HMI-ENT-1161  
14 cm/5 1/2"

BRUENINGS  
HMI-ENT-1162  
19 cm/7 1/2"

HMI-ENT-1163  
19 cm/7 1/2"

WIEDER  
14 cm/5 1/2" 38 mm HMI-ENT-1164  
15 cm/6" 45 mm HMI-ENT-1165

Beveled edge




OBWEGESER  
18 cm/7"




HMI-ENT-1166  
6 mm



HMI-ENT-1167  
7 mm



HMI-ENT-1168  
9 mm



HMI-ENT-1169  
11 mm

Flat edge



OBWEGESER  
18 cm/7"



HMI-ENT-1170  
6 mm



HMI-ENT-1171  
7 mm



HMI-ENT-1172  
9 mm



HMI-ENT-1173  
11 mm




BOIES  
HMI-ENT-1174  
19.5 cm/8"




Ramus periosteal stripper




OBWEGESER  
20.5 cm/8"



HMI-ENT-1175  
5 mm



HMI-ENT-1176  
7 mm




HMI-ENT-1177  
9 mm


Ramus periosteal stripper V-notched



OBWEGESER  
17.5 cm/7"




HMI-ENT-1178  
7 mm



HMI-ENT-1179  
9 mm

Mandibular stripper/inferior border raspatory




HMI-ENT-1180  
left


10 mm




OBWEGESER  
18 cm/7"



HMI-ENT-1181  
right



HMI-ENT-1182  
5 x 16 mm



HMI-ENT-1183  
7 x 25 mm




HMI-ENT-1184  
10 x 35 mm



HMI-ENT-1185  
10 x 42 mm




HMI-ENT-1186  
12 x 55 mm



HMI-ENT-1187  
14 x 70 mm



HMI-ENT-1188  
16 x 80 mm



HMI-ENT-1189  
5 x 16 mm



HMI-ENT-1190  
7 x 25 mm




HMI-ENT-1191  
10 x 35 mm



HMI-ENT-1192  
10 x 42 mm




HMI-ENT-1193  
12 x 55 mm



HMI-ENT-1194  
14 x 70 mm



HMI-ENT-1195  
16 x 80 mm





OBWEGESER  
22 cm/8 3/4"



OBWEGESER  
22 cm/8 3/4"


Mandibular channel retractors




Kinnhalter  
Chin retractor



Mandibular retractors




Ramus retractor




OBWEGESER  
24.5 cm/9 $\frac{3}{4}$ "

Notched alveolar retractor




OBWEGESER  
HMI-ENT-1204  
22 cm/9"

Nasal spine retractor



OBWEGESER  
HMI-ENT-1205  
24.5 cm/9 $\frac{3}{4}$ "


Osteotomes for alveolar process



3 mm HMI-ENT-1206  
4 mm HMI-ENT-1207  
6.5 mm HMI-ENT-1208

OBWEGESER  
15.5 cm/6"


Osteotomes



8 mm HMI-ENT-1200  
11 mm HMI-ENT-1210

OBWEGESER  
21 cm/8 1/4"


Maxillary mobilizer



HMI-ENT-1211  
12 mm

OBWEGESER  
23 cm/9"


convex / concave



HMI-ENT-1212  
12 x 30 mm

Bone contouring forceps


concave / convex



HMI-ENT-1213  
13.5 x 30 mm

Bone contouring forceps

flat / flat




HMI-ENT-1214  
13.5 x 30 mm

OBWEGESER  
17 cm/6 3/4"

OBWEGESER  
17.5 cm/7"

OBWEGESER  
17.5 cm/7"

Extra oral Signotch Retractors




for right cheek

for left cheek


for right cheek

for left cheek

26,5 mm x 16 mm



25,5 mm



for LeFort osteotomies



½


PETRI-PTERYGOID  
HMI-ENT-1223  
23 cm/9"



½

PETRI-PTERYGOID  
HMI-ENT-1224  
23 cm/9"

12,5 mm / 27 mm



for vertical osteotomies






½

LEAVASSEUR-MERRIL  
HMI-ENT-1225  
25 cm/9¾"




½

LEAVASSEUR-MERRIL  
HMI-ENT-1226  
25 cm/9¾"



|        |              |
|--------|--------------|
| Length | NEGUS        |
| 85 mm  | HMI-ENT-1229 |
| 140 mm | HMI-ENT-1230 |



HMI-ENT-1231  
15 cm/6"



HMI-ENT-1232  
21 cm/8 1/4"



HMI-ENT-1233  
26 cm/10 1/4"

Thank you for going through Catalog



[www.healermedinstuments.com](http://www.healermedinstuments.com)

Email: [Healermedins@gmail.com](mailto:Healermedins@gmail.com)

Whatsapp & Call: +92 3034515151